Product specification

Issue Date:

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*Material Assembly Drawings in appendix 1

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8. Biological Requirements	
Specification	Reference
Biocompatibility The sterile products shall fulfil the requirements according to ISO 10993-1 and following FDA Blue Book-Memorandum G95-1	ISO 10993-1 et FDA Blue book memorandum G95-1
Residuals after EtO sterilization	ISO 10993-7
Pyrogenicity: LAL method with a maximum of 0,5 EU/ml	FDA Guideline on Validation of the Limulus Amebocyte Lysate Test as an End-Product Endotoxin Test for Human and Animal Parenteral Drugs, Biological Products, And Medical Devices”, dated December 1987

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11. Regulatory Information	
Europe CE-marked Medical device class	Medical Directive 93/42/EEC Classe IIb

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

Stability test Apomorphine

To whom it may concern

The stability test has been carried out in order to evaluate possible changes in the stability of Apomorphin hydrochloride (Apo-go®) during normal use of the infusion sets neria™, neria™ detach and neria™ soft in a 12 hours exposure test.

The stability tests were carried out according to the experimental conditions described in the test report, documenting that the infusion sets neria™, neria™ detach and neria™ soft do not have any negative influence on the stability of Apomorphin hydrochloride (Apo-go®) during the test period for up to 12 hours.¹

Furthermore, a function test was conducted to detect any physical changes in neria™, neria™ detach and neria™ soft infusion sets after being exposed to the medication and to a test environment of 37 ° C. The function test was passed and no further test or actions are needed.¹

Christina Bach Pilegaard Søholm
Regulatory Affairs Coordinator
Unomedical a/s

Aaholmvej 1-3, Østed
Denmark

Christina B.P. Søholm
Regulatory affairs Coordinator
Unomedical a/s
Åholmvej 1-3, Østed
4320 Lejre · Denmark
Phone: +45 4642 7962
Fax: +45 4642 7815

8/7-2010

Date

neria is a trademark of Unomedical a/s
All other trademarks are the property of their respective owners

© 2009 Unomedical a/s
AP-008870-MM

Reference:

1 Stability test and function test of Comfort and Contact Detach by means of the medicament APO-go. Analysis report no. 294690. January 14-April 1, 2009. Data on file. Unomedical.

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